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- [1] PDB, Protein Data Bank: www.expasy.org, adiponectin: Q 15848 (21 oct.2005)
- [2] Meier U and Gressner AM (2004) Endocrine Regulation of Energy Metabolism: Review of Pathobiochemical and Clinical Aspects of Leptin, Ghrelin, Adiponectin, and Resistin. *Clin Chem* 50:1511-1525
- [3] Tsao TS, Lodish HF and Fruebis J (2002) ACRP30, a new hormone controlling fat and glucose metabolism. *Eur J Pharmacol* 440:213-221
- [4] Shapiro L and Scherer PE (1998) The crystal structure of a complement-1q family protein suggests an evolutionary link to tumor necrosis factor. *Curr Biol* 12:335-338
- [5] Bobbert T, Rochlitz H, Wegewitz U, Akpulat S, Mai K, Weickert MO, Mohlig M, Pfeiffer AF and Spranger J (2005) Changes of adiponectin oligomer composition by moderate weight reduction. *Diabetes* 54:2712-2719
- [6] Lara-Castro C, Luo N, Wallace P, Klein RL and Garvey T (2006) Adiponectin Multimeric Complexes and the Metabolic Syndrome Trait Cluster. *Diabetes* 55: 249-259
- [7] Duntas LH, Popovic V and Panotopoulos G (2004) Adiponectin: Novelities in Metabolism and Hormonal Regulation. *Nutr Neurosci* 7:195-200
- [8] Tsao TS, Lodish HF and Fruebis J (2002) ACRP30, a new hormone controlling fat and glucose metabolism. *J Biol Chem* 50:50810-50817
- [9] Nakashima R, Kamei N, Yamane K, Nakanishi S, Nakashima A and Kohno N (2006) Decreased Total and High Molecular Weight Adiponectin Are Independent Risk Factors for the Development of Type 2 Diabetes in Japanese-Americans. *Clin Endocrinol Metab* 91:3873-3877
- [10] Hara K, Horikoshi M, Yamauchi T, Yago H, Miyazaki O, Ebinuma H, Imai Y, Nagai R, and Kadowaki T (2006) Measurement of the High-Molecular Weight Form of Adiponectin in Plasma Is Useful for the Prediction of Insulin Resistance and Metabolic Syndrome. *Diabetes Care* 29:1357-1362
- [11] Nakano Y, Tobe T, Choi-Miura NH, Mazda T and Tomita M (1996) Isolation and Characterization of GBP28, a Novel Gelatin-Binding Protein Purified from Human Plasma. *J Biol Chem*

120:803-812

[12] Bodles AM, Banga A, Rasouli N, Ono F, Kern PA and Owens RJ (2006) Pioglitazone increases

secretion of high-molecular-weight adiponectin from adipocytes. *Am J Physiol Endocrinol Metab* 291:E1100-E1105

[13] Halperin F, Beckman JA, Patti ME, Trujillo ME, Garvin M, Creager MA, Scherer PE and Goldfine

AB (2005) The role of total and high-molecular-weight complex of adiponectin in vascular function in offspring whose parents both had type 2 diabetes. *Diabetologia* 48:2147-2154

[14] Fisher FF, Trujillo ME, Hanif W, Barnett AH, McTernan PG, Scherer PE, and Kumar S (2005)

Serum high molecular weight complex of adiponectin correlates better with glucose tolerance than total serum adiponectin in Indo-Asian males. *Diabetologia* 48:1084-1087

[15] Wang Y, Lu G, Wong WP, Vliegenthart JF, Gerwig GJ, Lam KS, Cooper GJ and Xu A (2004) Proteomic and functional characterization of endogenous adiponectin purified from fetal bovine serum. *Proteomics* 4:3933-3942

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